


other airway devices  
[created by Paul Young  
22/12/07]

larynx tube



- Sizes
- Size 3, for Child and Small Adult (Yellow Connector)
  - Size 4, for Medium Adult (Red Connector)
  - Size 5, for Large Adult (Violet Connector)

insertion:



combitube



features:

- white normal ETT
- Blue opening distal of the large cuff .
- Indicators for intubation depth.
- Cuff blue: 100ml air, seals the hypopharynx.
- Cuff white: 15 ml air, seals trachea / or esophagus
- can be inserted blindly

disadvantages:

- (i) Mostly adult sizes
- (ii) Potential for esophageal trauma
- (iii) Problems maintaining seal in some patients

general:

- Cuffed oropharyngeal airway
- Easy to use (like guedel)
- no protection against aspiration
- No protection against gastric inflation

COPA tube

